

clar N

WIRELESS HIGH PERFORMANCE LED HEADMIRROR

- Wireless and balanced headband with built in battery compartment
- Lifetime of the LED (140 lumen) up to 50.000 hours
- Significantly less energy consumption and heat generation
- > Built-in battery compartment and integral charge jack, optionally fitted with 2 lithium batteries size CR 123 A, 4 batteries size AAA (alkaline) or 4 rechargeable batteries size AAA (NiMH)
- > Battery operation approx. 4 hrs. with lithium batteries, approx. 90 min. with alkaline batteries Type AAA
- > Rechargeable battery operation: approx. 90 min. with fully charged battery
- > Headband cleaning made easy by inner, removable and washable padding
- > LED, 1 Watt with 4.000 - 4.200 Kelvin, CRI = 95 or with 6 V vacuum lamp
- > Weight:

without batteries	with CR 123 A lithium batteries	with rechargeable AAA batteries
209 g / 0.46 lb	242 g / 0.54 lb	258 g / 0.57 lb
- > Diameter mirror: 55 mm / 2.16 in



With 1 set of CR 123 A lithium batteries

6 V vacuum **No. 6070**
LED **No. 6072**

With NiMH rechargeable AAA batteries and plug-in charger

6 V vacuum 230 V **No. 6071**
6 V vacuum 120 V **No. 6074**
LED 230 V **No. 6073**
LED 120 V **No. 6075**

Replacement lamps

6 V vacuum 6 pcs. **No. 11301**
LED 1 item **No. 11302**

Batteries

Lithium CR 123A 2 pcs. **No. 11287**
AAA NiMH rechargeable batteries 4 pcs. **No. 11288**

Plug in chargers

230 V **No.:11284**
120 V **No. 11286**

Headband

No. 11295

Headlights only

6 V vacuum lamp **No. 11280**
LED Lamp **No. 11281**

Foam-rubber band

120 mm long **No. 11289**

Fabric points

Ø 35 mm (1.3 in) **No. 11293**
Ø 47 mm (1.8 in) **No. 11294**

Battery compartment cover

for lithium CR 123A **No. 12680**
for AAA NiMH rechargeable batteries **No. 12681**

Distance	400 mm / 15.75 in	200 mm / 7.87 in
6 V vacuum	4.000 lux	7.300 lux
LED	11.600 lux	30.000 lux

ziegler

HEAD MIRROR

- > Mirror Ø 90 mm (3.54 in)
- > Sturdy ball-and-socket joint - easy to assemble
- > Adjustable headband

Headlight **No. 6005**
Headband **No. 11305**
Headlights only **No. 11312**



901004 Rev. N 2023-03-27

Please visit riester.de for more information